SINCE 1966
MADE IN ITALY

Smarty

by **FLAEM**®

MANUALE
ISTRUZIONI
D'USO

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INSTRUCTIONS
FOR USE
MANUAL

ENGLISH
Page 19

MODE
D'EMPLOI

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NEDERLANDS
Page 55

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
MANUAL DE
INSTRUCCIO-
NES DE USO

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Assembly diagram - Schéma de montage -
Verbindungsschema - Anschlusschema -Cro-
quis de conexiones- Schemat połączeń dla
nomenklatury

Smarty

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



Portable device

for aerosol therapy using

the latest in PMVT technology (Passive MESH Vibrating Technology)

Thank you for your purchase. Our goal is the complete satisfaction of our customers by offering cutting-edge products for the prevention and treatment of the respiratory tract ailments.


Visit our website www.flaem.it to view the whole range of Flaem's products. Read these instructions carefully in order to use the device correctly.

We recommend that you store this manual for future consultation.

Smarty is a medical device for home use to nebulise and administer medication


Intended use: this is an electronic Medical Device for aerosol therapy. **Smarty** is suitable for nebulising **medications which are either in solution** eg mucolytics, or **suspension** eg corticosteroids such as **beclomethasone dipropionate** and **budesonide**) and saline solution, prescribed or recommended by your doctor upon assessing the patient's general conditions. This unit may be used in medical facilities (hospitals, clinics etc.) and at home.

WARNING
TO PREVENT DAMAGE TO THE MESH,
CAREFULLY FOLLOW STEPS 1 TO 4 AS LISTED BELOW


For more information, a video tutorial is available at smarty.flaem.it

DEVICE COMPONENTS


1- Main Unit


2- Removable medication chamber

2G- Nebuliser head


6 - 2 1.5V AA alkaline batteries


3- Mouthpiece

4- Child mask

5- Adult mask

1A - ON DEMAND nebulisation button

1B - CONTINUOUS nebulisation
START AND STOP button

1C - Medication chamber release button

2A - Medicine cup


2B - Nebulising chamber

2C - Medication cup locking clip

2D - Medication cup cover

2E - Nebulising chamber and
Medication cup release clips

2F - Accessories connection port


IMPORTANT WARNINGS

- As with any electrical device, especially in the presence of children, the use of **Smarty** requires some basic safety precautions.
- The device should not be used as a toy. Please take particular care when it is used by children.
- If the device is used by or in the presence of children or people who require special assistance, the supervision of an adult who has read this manual is necessary.
- This product must not be used by unconscious patients or not breathing spontaneously
- Some parts of the device are small enough to be swallowed by children; keep the device out of reach of children.
- Never use batteries other than those recommended in this manual.
- Use only batteries of known brands. Always replace both with new ones and do not mix new and old batteries. AA alkaline or Lithium batteries and rechargeable Ni-Mh can be used. Rechargeable batteries shall be charged using a battery charging device that is not supplied with this device.
- The average battery life depends on the brand used.
- It is very important to use batteries of the same brand and same type.
- Remove the batteries before storing the unit when it is not used for long periods.
- Smarty should not be used in the presence of an anaesthetic mixture which is flammable with oxygen or nitrous oxide.
- **The device housing has been designed for limited contact with liquids (no immersion). Do not wash the device under running water or by immersion. Carefully follow the "CLEANING, SANITISATION, DISINFECTION" in this manual. Absolutely do not use alcohol to clean the device**
- Do not use the device while taking a shower or a bath.
- Do not operate the device while driving or in any other situation where a distraction could be dangerous for the user, for close people or animals or surrounding objects.
- Do not expose the device and the batteries to extreme temperatures. Keep the device and batteries away from heat sources, direct sunlight or excessively hot and humid environment.
- The average duration for the device used for about 20 min a day with medicinal products (4 applications of 5 min) is set out below:
 - Main unit 5 years
 - Masks and Mouthpiece 5 years
 - Medication chamber 1 year
 - Nebuliser head 1 year (Component subject to wear)
- Should your device fail to provide the expected performance, contact the authorised service centre for clarifications.
- Repairs must be performed by authorised FLAEM personnel only, by complying with the information provided by the manufacturer. Any

unauthorized repairs will void the warranty and may pose a safety hazard for the user.


-  **WARNING** Do not modify this device without authorisation of the manufacturer.
- The Manufacturer, Importer and Seller will be held responsible for the safety, reliability and performance, only if the device is used in accordance with the instruction manual.
- Interactions: the materials used have been tested according to the standards of biocompatibility (ISO 10993-5 and ISO 10993-10) in compliance with the essential requirements of the 93/42 EEC Medical Device Directive. The materials used in the device are biocompatible in accordance with the provisions of Directive 93/42 EC and subsequent amendments. However, the possibility of occurrence of allergic reactions cannot be entirely excluded.
- If the device is used with different medications, it is recommended to remove all residues. Therefore, clean it after each treatment so as to obtain the highest degree of hygiene and to optimise the life and operation of the device
- The device does not contain parts to be repaired by the user. The warranty does not cover batteries or damage caused by incorrect, used up or improperly stored batteries.
- The manufacturer should be contacted for reporting problems and/or unexpected events related to device operation, and for any clarifications on use, maintenance/cleaning.
- Only use original accessories and spare parts by Flaem, we disclaim any liability in the event of using non original spare parts or accessories

OPERATING INSTRUCTIONS

1 - BATTERY INSTALLATION

- 1.a Use a small coin to open the **Battery Compartment cover (1D)** located on the bottom of the unit.

1.b Insert the two **AA batteries (6)**, in the **Battery Compartment (1E)**, making sure that the polarity is correct.


1.c Close the **Battery Compartment Cover (1D)** making sure that the cover is properly.

2 - PREPARATION

WARNING!

Before each use, the unit and the accessories should be cleaned and/or disinfected according to the **DISASSEMBLY, CLEANING, SANITISATION, DISINFECTION** instructions (paragraph 5/6)

We recommend personal use, both of the nebuliser and accessories to avoid possible risk of infection. **Use only FLAEM original accessories.**

Wash your hands thoroughly before preparing the medication

2.1 Lift the **locking clip (2C)** of the **Medication cup cover (2D)**

2.2 Open the **Medication cup cover (2D)** completely

2.3 Pour the medicine (MAX 8 ml) into the **Medication cup (2A)**

2.4 Close the **Medication cup cover (2D)**

2.5 Fasten the **locking clip (2C)**

WARNING!


After pouring the medicine into the **Medication cup (2A)**, the medicine must be nebulised.

3 - UNIT CONFIGURATION

3.1 Connect the desired accessory to the Accessories **connection port (2F)**:


- **Mouthpiece (3)**
- **Child mask (4)**
- **Adult mask (5)**

If a mask is used, place it on the face as shown in the picture

4 - OPERATION

For effective therapy, after setting the appliance up, sit in a comfortable and relaxed position and proceed. Select the preferred delivery mode ("**ON DEMAND**" or "**CONTINUOUS**")

4.1 By pressing buttons **(1A-1B)** they light up as follow:


FLASHING BLUE: for a few seconds, until a complete initialization of device.

STEADY BLUE: nebulisation

FLASHING RED: batteries are used up; replace with new batteries of the same type

4.2 ON DEMAND button (1A) Press this button when inhaling to obtain the delivery of medication without waste. The device only works if you hold down the button, which will turn blue.

4.3 START AND STOP / CONTINUOUS button (1B): Press this button to obtain the "**CONTINUOUS**" delivery of medication (fast treatment).

To switch off the unit, press the **button (1B)** again.


4.4 Sit comfortably holding the device in your hand, place the **Mouthpiece (3)** in your mouth then breath in and out slowly through the mouth. Alternatively, you can use the **Child Mask (4)** or the **Adult Mask (5)**.

4.5 To increase the effectiveness of treatment, breathe slowly and deeply and after inhaling hold your breath for a moment, so that the inhaled medicine can deposit along the respiratory tract. Then exhale slowly.


4.6 The unit turns off automatically:

4.6.a after 5 minutes of use

4.6.b when, during use, the **medication chamber (2)** is accidentally removed .

4.7 If you do not use the appliance for a pause longer than 5 minutes, place it down on a surface as shown in the picture. At the end of each use store the device complete with accessories in a dry place away from dust.


When treatment is finished proceed as indicated in paragraphs 5 and 6.


5 - DISASSEMBLY

When treatment is finished:

- 5.1 Press the **Medication chamber release button (1C)** positioned on the back of the unit.
- 5.2 Remove the **Medication chamber (2)**.




- 5.3 To access the **Nebuliser head (2G)** in the **Nebulising chamber (2B)** release the **Medication cup locking clip** and open the **Medication cup cover (2D)**.
- 5.4 Press simultaneously the two **release clips (2E)** on both sides of the Medication chamber.
- 5.5 Rotate until the **Nebulising chamber (2B)** opens completely.

- 5.6. Remove the nebuliser head (2G), while at the same time pressing the side fins.

⚠ WARNING

- Do not use tools: this may cause irreparable damage to the **Nebuliser head (2G)**.
- Do not press the centre of the **Nebuliser head (2G)** with your fingers or tools from the inside: this may cause irreparable damage.


5.7 Carefully clean it according to the **CLEANING, SANITISATION, DISINFECTION** instructions in paragraph 6.

6 - CLEANING, SANITISATION, DISINFECTION


Before and after each use, the device and the accessories must be cleaned and disinfected properly as described below. If this is not done, some microorganisms may collect in the device, causing the risk of infection. Do not use alcohol or other solvents to clean the device.

6.1. CLEANING Accessories and Removable Upper Part (2):


6.1.1. WASHING: Wash the components **2-2G-3-4-5** with lukewarm drinking water with dish washing soap (non-abrasive), then rinse thoroughly with a jet of lukewarm drinking water to remove any detergent residue.


6.1.2. DRYING

6.1.2.a. After washing and rinsing all the components, shake them to remove excess water and reassemble them, as described in the following paragraph **7. REASSEMBLY OF THE DEVICE**. At this point, to eliminate the residual water that has deposited in the **nebuliser head (2G)**, press **Button 1B** (side figure) and wait for the device stop to nebulizing. Turn off the unit. This operation is very important, because it prevents the deposit of limestone in the micropores of the **Nebuliser Head (2G)** which could affect proper nebulisation of the drug in subsequent therapeutic applications.


6.1.2.b. To complete the drying of all the components, leave them out and in a dry place (e.g. not in the bathroom) as into the side figure.

WARNING

Do not clean the **Mesh (2G)** using brushes or other instruments: this may cause irreparable damage.

6.2 CLEANING the Main Unit


To clean the **Main Unit (1)** use a cloth slightly moistened with an antibacterial detergent, then dry it with a soft paper towel.

WARNING

Do not immerse the **Main Unit (1)** in water and do not wash it under running water.

6.3. SANITISATION






6.3.1. Immerse for at least 30 minutes the opened **Removable Medication Chamber (2)** and the accessories **2G-3-4-5**, in a solution of 50% water and 50% white vinegar, then rinse thoroughly with hot running water.

6.1.2. DRYING

6.1.2.a. After washing and rinsing all the components, shake them to remove excess water and reassemble them, as described in the following paragraph **7. REASSEMBLY OF THE DEVICE**. At this point, to eliminate the residual water that has deposited in the **nebuliser head (2G)**, press **Button 1B** (side figure) and wait for the device stop to nebulizing. Turn off the unit. This operation is very important, because it prevents the deposit of limestone in the micropores of the **Nebuliser Head (2G)** which could affect proper nebulisation of the drug in subsequent therapeutic applications.

6.1.2.b. To complete the drying of all the components, leave them out and in a dry place (e.g. not in the bathroom) as into the side figure.

DO NOT IMMERSE IN
SOLUTIONS BASED ON
CHLORINE

6.4.1 METHOD A: DISINFECTION BY STEAM

6.4.1.a Disinfect the accessories provided in Table 1 following the disinfection procedure with device for feeding-bottle (without microwave), it requires a treatment duration of at least 15 minutes. It is necessary to use demineralised or distilled water to prevent limestone deposits, which could impair the function of the nebuliser head. Carefully read the instructions for use of the device in question, paying special attention to the amount of needed water and the method to be followed to carry out the disinfection phases.

6.4.2 METHOD B: DISINFECTION BY BOILING

6.4.2.a Disinfect the accessories provided in Table 1 by boiling for 20 minutes; it is necessary to use de mineralized or distilled water to prevent limestone deposits, which could impair the function of the nebuliser head.

6.4.2.b DRYING: Before making a new application wait until the components have cooled down and dried.

6.4.3 METHOD C: DISINFECTION BY CHEMICAL


WARNING Except the nebuliser head
2G

6.4.3.a The disinfectant to be used must be an electrolytic chloroxidizer (active principle: sodium hypochlorite), specific for disinfecting, which is available in any pharmacy.

6.4.3.b Fill a container, whose size is suitable to hold all of the accessories to be disinfected listed in Table 1, with a drinking-water and disinfectant solution, according to the proportions indicated on the packaging of the disinfectant itself.

6.4.3.c Completely immerse each accessory in the solution, taking care to avoid the formation of air bubbles that may come in contact with the accessories. Leave the accessories immersed for the amount of time indicated on the packaging of the disinfectant, and associated with the concentration chosen to prepare the solution.

6.4.3.d Remove the accessories and rinse


them with plenty of lukewarm drinkable water.

6.4.3.e DRYING: carefully perform all the operations described in paragraph 6.1.2. of the user manual.

6.4.3.f Dispose of the solution following the instructions provided by the disinfectant manufacturer.

METHODS DISINFECTABLE PARTS		method a) Equipment bottle disinfection by steam (not a microwave oven)	method b) Boiling with demineralised or distilled water	method c) Electrolytic chloroxidizer	method d) Alcohol
2G Nebulizer head		YES	YES	NO	NO
3 Mouth piece		YES	YES	YES	YES
4 –5 Adult and Child mask		YES	YES	YES	YES
2 Complete medication chamber (2G excluding)		YES	YES	YES	NO


7 - DEVICE RE-ASSEMBLY


7.1. Insert the Spraying Head (2G) in its seat on the front of the medicinal product container until the double click sound indicates successful coupling.

WARNING


Do not use any tools, you could damage the head.


7.3. Close the complete pull-out medicine container (2) by rotating the Front Part (2B) towards the Drug Container (2A) until the two Container Buttons (2E) are hooked.

7.4. Close the **Medication Cup Cover (2D)** by means of the **locking clip (2C)**

7.5. Insert the **Removable Medication container (2)** onto the **Main Unit (1)** matching the guides on the two parts.

7.6. Push the **Removable Medication container (2)** right down until the **Medication container release button (1C)** located in the rear of the **Main Unit (1)** springs with a **CLICK**.

TROUBLE-SHOOTING

PROBLEM	CAUSE	SOLUTION
The device switches on but does not nebulise or nebulises little.	Malfunction of the Mesh (2G)	Soak the Nebulising head (2G) in a solution of water (50%) and white wine vinegar (50%), rinse and fit back on as described in paragraph 6.3 SANITISING
	The nebuliser head (2G) is broken	Replace it with a new Nebuliser head (2G)
	No medicine supply in the nebulising area	Verify that the small hole in the centre of the Medication cup cover (2D) is not obstructed. If so clean it with a needle
	The nebulising speed depends on the type of medicinal product used	The inhaling treatment time depends on the type of medicinal product and on the patient's inhaling ability
The device switches off after a few seconds (blue light).	No electrical contact between the Removable Medication chamber (2) and the Main Unit (1)	Check that the Removable Medication chamber (2) is properly hooked to Main Unit (1)
The device switches off	It's been 5 minutes since device switching on	Turn on the unit again
The START AND STOP Button (1B) shows a flashing red light and the device does not nebulise or switches off.	Batteries are used up or low	- Replace the batteries if they are alkaline. - Recharge the Batteries if they are rechargeable Ni-MH or Lithium type.
The device does not switch on.	No electric power supply	- Verify that Batteries are inserted with the correct polarity. - Used up or low Batteries. Replace or recharge.

If after checking the aforementioned conditions, the unit still does not nebulise, we recommend that you contact the dealer or an authorised Flaem service centre.

8 - TECHNICAL CHARACTERISTICS

Smarty Electronic Nebuliser Model P0315EM

Power Supply	3.0Vdc with 2 Batteries AA Alkaline, 1.5Vdc 2.4Vdc with 2 Batteries AA Rechargeable Ni-MH or Lithium, 1.2Vdc
Power	1.5W
Transducer frequency	170 KHz
Noise level (at 1 m)	less than 20 dBA
Device dimensions	54(W) X 69(D) X 140(H) mm
Device weight (excluding batteries)	125g
Carry bag dimensions	140(W) X 80(D) X 150(H) mm
Delivery rate ml/min⁽¹⁾	≥ 0.25 ml
Characterization * MMAD	4,00 µm approx
Breathable fractions < 5 µm (FPF):	65% approx
Medication cup capacity⁽¹⁾	8 ml

⁽¹⁾The nebulising speed has been measured with a 0.9% saline solution at 23° C according to the internal Flaem procedure I29-P07.5. It may vary according to the nebulising head supplied, medicinal product and environmental conditions. Values shown for the aerosol delivery may also vary depending on the patients respiratory Capacity.

* Please visit Smarty.Flaem.it or contact Flaem Nuova S.p.A Via Colli Storici, 221 25015 S. Martino della Battaglia (BS) Italy - Tel +39 030 9910168 for information concerning nebulisation characteristics.

APPLIED PARTS

The BF type applied parts are: patient accessories (3,4,5)

9 - OPERATING CONDITIONS

Temperature	min 10°C ~ max 40°C
Air humidity	min 10% ~ max 75% RH
Atmospheric Pressure	min. 690 hPa ~ max. 1060 hPa

10 - STORAGE CONDITIONS

Temperature	min 10°C ~ max 35°C
Air humidity	min 10% ~ max 75% RH
Atmospheric Pressure	min. 690 hPa ~ max. 1060 hPa

If the appliance is stored at a temperature other than indicated under the storage conditions, leave the device at room temperature for at least 1 hour prior to use. With medicinal products in suspension or particularly viscous, the information provided according to standard EN13544-1 might be subject to changes.

Temperature, atmospheric pressure and humidity might affect the appliance's performance.

11 - SYMBOLS

	Warning			Battery disposal: Used up batteries must be disposed of in appropriate containers.
	Type BF applied parts			
	WARNING! Check the instructions for use manual			Manufacturer
	Direct current			Keep dry
				Serial number of the device
	CE Marking medical ref. Dir 93/42 EEC Directive and subsequent updates			Positive batteries polarity
				Negative batteries polarity

12 - DISPOSAL OF DEVICE

 In conformity with Directive 2012/19/EC, the symbol shown on the device to be disposed of indicates that it is considered as waste and is therefore subject to "sorted waste collection". The user must therefore take (or have taken) the above waste to a pre-sorted waste collection centre set up by the local authorities, or else give it back to the dealer when purchasing a new appliance of the same type. Pre-sorted waste collection and the subsequent treatment, recovery and disposal operations favor the production of appliances made of recycled materials and limit the negative effects of any incorrect waste management on the environment and public health. The unlawful disposal of the product by the user could result in administrative fines as provided by the laws transposing Directive 2012/19/EC of the European member state or of the country in which the product is disposed of.

13 - ELECTROMAGNETIC COMPATIBILITY

This device has been designed to satisfy requirements currently required for electromagnetic compatibility (EN 60 601-1-2:2007). Electrical medical devices require special care during installation and use with respect to EMC requirements. It is therefore required that they be installed and/or used according to the manufacturer's specification. Potential risk of electromagnetic interference with other devices. Radio and mobile telecommunications devices or portable RF (mobile phones or wireless connections) may interfere with the operation of electrical medical devices. For further information visit www.flaemnuova.it website. The Device may be subject to electromagnetic interference if other devices are used for specific diagnosis or treatments. Flaem reserves the right to make technical and functional modifications to the product without notice.

14 - EXPECTED DURATION

Expected duration: The service life of the various components set out below refer to using the device for about 20 min a day (4 applications of 5 min.)

Main unit	5 years
Masks and Mouthpiece	5 years
Medication chamber	1 year
Nebuliser head	1 year
(Component subject to wear)	

AVAILABLE SPARE PARTS AND ACCESSORIES


Mouthpiece
(made of silicone)

ACO563P

Child mask
(Soft-touch silver)

ACO438P


Adult mask
(Soft-touch silver)

ACO437P

Nebuliser head

ACO655P

**Complete removable
medication container**

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© 2015 FLAEM NUOVA® All rights reserved
FLAEM NUOVA S.p.A.